



5°C Pharmacy Refrigerator STPR-A10

Overview

5°C Pharmacy Refrigerator STPR-A10 has a capacity of 226L and maintains a temperature range of 2°C to 8°C. It features a suction circulation system, ensuring uniform cooling and consistent temperature distribution. Our pharmacy refrigerator uses a CFC-free refrigerant and an energy-saving system, ensuring eco-friendly and efficient operation. It is equipped with universal castor wheels for easy mobility and stable positioning.

Features :

Forced air cooling

Lockable doors

Fast pull down

Audible and visual alarms

Variable space

Specifications :

Capacity	226 L
Temperature Range	2°C to 8°C
Ambient Temperature	10°C to 32°C
Cooling Method	Forced air cooling
Defrost Mode	Auto
Refrigerant	R600a
Sensor	NTC
Compressor	1
Temperature Controller	Microprocessor
Display	Digital display
External Material	SPCC sprayed steel
Internal Material	HIPS
Alarms	High/low temperature Power failure Sensor error Controller malfunction Door ajar 8h Power failure backup
Door type	Heated glass door
Shelf	4
Noise	55db
Power	221W
Current	1.4A
Power Consumption	3.47 kWh/24h
Power Supply	220 V 50 Hz/110V 60Hz
Internal Dimension	534 × 406 × 1150 mm
External Dimension	625 × 574 × 1707 mm
Packaging Dimension	707 × 676 × 1745 mm
Net Weight	73 Kg
Gross Weight	76 Kg

Standard Accessories :

1/25mm access port

Castors

Optional Accessories :

Data logger

USB

Remote Alarm port

Applications :

5°C Pharmacy Refrigerator STPR-A10 is commonly used in pharmacies, hospitals, and clinics to store temperature-sensitive medications and vaccines safely.

Standard Accessories :

1/25mm access port

Castors

Optional Accessories :

Data logger

USB

Remote Alarm port

