



5°C Pharmacy Refrigerator STPR-C14

Overview

5°C Pharmacy Refrigerator STPR-C14 provides 130L capacity and maintains temperatures between 2°C and 8°C. It uses forced air cooling and a high-performance compressor for stable temperature control. Safety features of our pharmacy refrigerator include temperature alarms, sensor failure alerts, power outage warnings, and a door open alarm. The heated double-layer glass door improves insulation and visibility.

Features :

Lockable doors

Universal wheels

Audible and visual alarms

Variable space

LED lighting

EBM silent fan

CFC-free refrigerant

Specifications :

Capacity	130 L
Temperature Range	2°C to 8°C
Refrigeration mode	Forced air cooling
Defrost Mode	Auto
Refrigerant	R134a
Foaming thickness	40 mm
Compressor	1
Temperature Controller	Microprocessor
Display	Digital display
Noise	52 dB
Doors	1
Shelf	4
Power	240 W
Power Consumption	2.8 kWh/24h
Power Supply	220 V/50 Hz
Net Dimension	432 × 500 × 843 mm
External Dimension	522 × 595 × 1365 mm
Packaging Dimension	610 × 690 × 1515 mm
Net Weight	62 Kg

Optional Accessories :

USB

Applications :

5°C Pharmacy Refrigerator STPR-C14 is commonly used in pharmacies, hospitals, and clinics to store temperature-sensitive medications and vaccines safely.

Optional Accessories :

USB